

PTR-MS Technology

Advanced Technology: Proton Transfer Reaction Mass Spectrometer (PTR-MS)


Mass spectrometry is a critical component of many analytical techniques used for the identification and quantification of organic compounds in a wide variety of sample matrices from both environmental and industrial sources. The foundational principle is that individual chemical compounds have a known molecular weight that can be used to identify and quantify them. The physics behind the analytical methods are beyond the scope of this article, but it is typically necessary to add a charge (ionize) to each molecule to allow for mass separation, detection and quantification using a mass spectrometer. Analysis using a mass spectrometer, commonly results in fragmentation of larger molecules into smaller ones, complicating the interpretation, as the pattern of the fragmented masses needs to be identified and used towards the identification and quantification of the original larger molecules. The delicate nature of these analytical operations requires that most mass spectrometers can only operate within the confines of a laboratory, analyzing discrete samples introduced into these instruments.

The advent of PTR-MS technology has addressed several of the limitations of laboratory benchtop units described above. PTR-MS is a “soft” ionization mass spectrometer. Charge is added to each molecule using water, oxygen, or nitrogen oxide molecules which suppress the fragmentation of the compounds, allowing for a more robust analysis of the mass spectrum. PTR-MS is particularly adept at identification and measurement of volatile organic compounds (VOCs) in air. It can provide real-time, continuous measurement of multiple volatile organic compounds (VOCs) at ultra-trace levels, i.e., in the parts per trillion (ppt) and low parts per billion (ppb) range. A truck-mounted solution also gives the instrument mobility for environmental monitoring and sourcing investigations.

The PTR-MS has been used in many peer-reviewed studies for the detection of VOCs. Studies include ambient measurements for the investigation of the air quality of urban and suburban areas.^{1,2} The information that the PTR-MS provides can be used in air quality source apportionment studies in order to gain insight of the sources of different air toxics such as toluene, xylene, benzene; biogenic compounds emitted by vegetation, such as isoprene, monoterpenes and sesquiterpenes; and other compounds such as acetaldehyde, formaldehyde and acetone.^{3,4} Other studies have used the PTR-MS for the detection of the emissions from different sources, such as biomass burning^{5,6}, cooking activities⁷ and diesel exhaust⁸. Other applications include the testing of emissions from materials, investigation of indoor air quality and vapor intrusion studies⁹⁻¹². Most recently the device was used for the breath analysis and identification of biomarkers of Covid-19 patients¹³. The ability of real time detection of VOCs makes the PTR-MS the ideal instrument for these applications.

Case Study 1: Laboratory testing - headspace sampling

Due to its high sensitivity and low time resolution, the PTR-MS is the perfect tool for headspace sampling. RJ Lee Group has implemented the PTR-MS in several headspace projects for the detection of emissions from carbon black, canned food products, as well as products such as gloves and face masks. In a previous study, RJ Lee Group investigated the compounds emitted by a medical device. The device was placed in a testing Tedlar bag, which was filled with clean air. The PTRMS was used to monitor the VOC concentration inside the bag. After the device was turned on, the concentration of 6 compounds increased. After two hours of testing, compounds such as butadiene and hexynes had reached a plateau at a concentration lower than 5 ppb (Figure 1). The concentration of acetone, acetic acid, isoprene and butene continued increasing during the whole testing period.



Testing of emission collected with Tedlar bags using the PTR-MS.

CONNECT WITH AN EXPERT

800.860.1775 | WWW.RJLEEGROUP.COM

 **RJ LEE GROUP**
DELIVERING SCIENTIFIC RESOLUTION




Figure 1:
Headspace sampling of the emissions generated by a medical device using a PTR-MS.

Case Study 2: Field measurements mobile applications

The PTR-MS can also be used to study the air quality of an area. The combination of the PTRMS with RJ Lee Group's mobile laboratory can allow for temporal and spatial resolution air quality measurements. Our scientists have conducted numerous fence-line monitoring and community air quality monitoring projects. In a recent study we measured the ambient concentration of benzene in the greater area of Philadelphia. We were able to identify the higher emitters and investigate their impact in the neighbor communities (Figures 2 and 3).





Figure 2:
The timeseries of benzene concentration during mobile measurements.



Benzene (ppb)

- 11.1
- 9.8
- 8.6
- 7.4
- 6.2
- 5.0
- 3.8
- 2.6
- 1.4
- 0.2

Figure 3:
The concentration of benzene in the greater Philadelphia area. The altitude and the color of the series is proportional to the concentration of benzene.

Applicable RJLG Services

- » Emission testing of materials and commercial products using the TGA-PTR-MS
- » Fenceline monitoring, industrial emissions evaluations
- » Environmental test chamber studies
- » Indoor air quality testing (IAQ)
- » Vapor intrusion studies
- » Source identification
- » Product off-gassing
- » Infection detection



The RJ Lee Group PTR-MS mobile laboratory at a Philadelphia area refinery.

References

1. De Gouw, J. A., Welsh-Bon, D., Warneke, C., Kuster, W. C., Alexander, L., Baker, A. K., Beyersdorf, A. J., Blake, D. R., Canagaratna, M., Celada, A. T., et al.: Emission and chemistry of organic carbon in the gas and aerosol phase at a sub-urban site near Mexico City in March 2006 during the MILAGRO study. *Atmos. Chem. Phys.*, 9, 3425-3442, 2009.
2. Shah, R. U., Coggon, M. M., Gkatzelis, G. I., McDonald, B. C., Tasoglou, A., Heinz Huber, Jessica Gilman, J., Warneke, C., Robinson, A. L. Presto, A.: Observational evidence of the importance of secondary organic aerosol of volatile chemical products, *Environ. Sci. Technol.*, 54 (2), 714-725, 2020.
3. Kaltonoudis, C., Kostenidou, E., Florou, K., Psichoudaki, M., Pandis, S. N.: Temporal variability and sources of VOCs in urban areas of the eastern Mediterranean, *Atmos. Chem. Phys.*, 16, 14825-14842, 2016.
4. Crippa, M., Canonaco, F., Slowik, J. G., El Haddad, I., DeCarlo, P. F., Mohr, C., Heringa, M. F., Chirico, R., Marchand, N., Temime-Roussel, B., Abidi, E., Poulain, L., Wiedensohler, A., Baltensperger, U., and Prévôt, A. S. H.: Primary and secondary organic aerosol origin by combined gas-particle phase source apportionment, *Atmos. Chem. Phys.*, 13, 8411-8426, 2013.
5. Tasoglou, A., Saliba, G., Subramanian, R., Pandis, S. N.: Absorption of Chemically Aged Biomass Burning Carbonaceous Aerosol, *J. Aerosol. Sci.*, 113, 141-52, 2017.
6. Gilman, J. B., Lerner, B. M., Kuster, W. C., Goldan, P. D., Warneke, C., Veres, P. R., Roberts, J. M., de Gouw, J. A., Burling, I. R., and Yokelson, R. J.: Biomass burning emissions and potential air quality impacts of volatile organic compounds and other trace gases from fuels common in the US, *Atmos. Chem. Phys.*, 15, 13915-13938, 2015.
7. Klein, F., Platt, S. M., Farren, N. J., Detournay, A., Bruns, E. A., Bozzetti, C., Daellenbach, K. R., Kilic, D., Kumar, N. K., Pieber, S. M., et al.: Characterization of gas-phase organics using proton transfer reaction time-of-flight mass spectrometry: cooking emissions, *Environ. Sci. Technol.*, 50, 1243-1250, 2016.
8. Inomata, S., Tanimoto, H., Fujitani, Y., Sekimoto, K., Sato, K., Fushimi, A., Yamada, H., Hori, S., Kumazawa, Y., Shimono, A., et al.: On-line measurements of gaseous nitro-organic compounds in diesel vehicle exhaust by proton-transfer-reaction mass spectrometry, *Atmos. Environ.*, 73, 195-203, 2013.
9. Sears, J.; Rogers, T.; McCoskey, J.; Lockrem, L.; Watts, H.; Pingel, L.; Conca, J. Proton Transfer Reaction Mass Spectrometry as a Real-Time Method for Continuous Soil Organic Vapor Detection. In *Continuous Soil Gas Measurements: Worst Case Risk Parameters*; ASTM International: 100 Barr Harbor Drive, PO Box C700, West Conshohocken, PA 19428-2959, 32-44, 2013.
10. Peinado, I., Marco, M., Franco, B., Scampicchio, M.: Hyphenation of Proton Transfer Reaction Mass Spectrometry with Thermal Analysis (TG/TR-MS) for Monitoring the Thermal Degradation of Retinyl Acetate, *Rapid Commun. Mass Spectrom.* 32, 57-62, 2018.
11. Rosales, C. M. F., Jiang, J., Lahib, A., Bottorff, B. P., Reidy, E. K., Kumar, V., Tasoglou, A., Huber, H., Dusanter, S., Tomas, A., Boor, B. E., Stevens, P. S.: Chemistry and human exposure implications of secondary organic aerosol production from indoor terpene ozonolysis, *Sci. Adv.*, 8, 2022.
12. Wu, T., Tasoglou, A., Huber, H., Stevens, P. S., Boor, B. E.: Influence of Mechanical Ventilation Systems and Human Occupancy on Time-Resolved Source Rates of Volatile Skin Oil Ozonolysis Products in a LEED-Certified Office Building, *Environ. Sci. Technol. Lett.*, 24, 16477-16488, 2021.
13. Liangou, A., Tasoglou, A., Heinz J. H. J., Wistrom, C., Brody, K., Menon, P. G., Bebekoski, T., Menschel, K., Davidson-Fiedler, M., DeMarco, K., Salphale, H., Wistrom, J., Wistrom, S., Lee, R. J.: A novel method for the identification of COVID-19 biomarkers in human breath using the PTR-ToF-MS, *EclinicalMedicine*, 42, 101207, 2021.

This summary description is offered for educational purposes only and no warranty (express or implied) is given on the applicability of this information to any particular part, material, application, or product. This article is not intended to be a complete dissertation on this topic and the reader is encouraged to contact RJ Lee Group for consultation regarding individual applications, materials, or situations.



Author: Antonios Tasoglou Ph.D., Senior Scientist

Dr. Antonios Tasoglou is a senior scientist in RJ Lee Group. He has 12 years of experience in air quality measurements of hazardous air pollutants - air toxics and particulate matter. He has conducted numerous field studies for the private and public sector. He has a deep understanding on the use of research grade, FEM and FRM instrumentation for the analysis of the air quality. He is an expert at PTR-MS applications and completed his doctoral studies at Carnegie Mellon University.

CONNECT WITH AN EXPERT
800.860.1775 | WWW.RJLEEGROUP.COM

 **RJ LEE GROUP**
DELIVERING SCIENTIFIC RESOLUTION